


Treadmill ergosprint

Treadmill



The figure could show not selected options and equipment ranges.



Treadmill ergosprint

Treadmill



The figure could show not selected options and equipment ranges.



Treadmill ergosprint

Technical Specifications

Device name	Treadmill ergosprint	Treadmill ergosprint P	
Control	via integrated interface no display, no keyboard	via user terminal with keyboard and display, integrated interface	
Running surface	Length: 59 in, Width: 20 in Access height: 9.06 in Running belt with slip resistant surface		
Speed range	0.3 – 11.0 mph (in steps of 0.06 mph) available at extra charge: 0 – 13.6 mph		
Acceleration	7 acceleration/deceleration levels between 131 s and 3 s from 0 to max. or from max. to 0		
Elevation	0 – 20% available at extra charge: 0 – 25%		
Motor system	2.2 kW (3 HP) 3-phase AC motor, maintenance free and brushless		
Power transmission	frequency inverter, poly-V-belt, very quiet operation		
Safety systems	CE0123; medical device directive 93/42/EEC, machinery directive 2006/42/EC; ISO 20957-1; EN 957-6; EN 14971; EN ISO 13485; IEC60601-1; EN 60601-1-2 (EMC approved); IEC 62304 emergency-off safety stop switch (mushroom push button for drive system power-off); emergency stop switch (safety lanyard with actuator, pull cord and clip); potential equalization bolt; transformer for potential-isolation from the mains		
Degree of protection	appliance class I/Type B/IP 00		
Classification	medical device risk class IIb according to MDD 93/42/EEC, active therapeutic medical device and active diagnostic medical device		
Digital interface	1 x RS 232 Com1 with 9600 bps. incl. PC-protocol h/p/cosmos coscom® v3/v4		
Voltage supply	200 – 240 Volt AC 1~/N/PE 50/60 Hz 16A, type C fuse; dedicated circuit line and protection		
Permissible load	max. 550 lbs		
Options			
Safety	Safety arch with chest belt and fall-stop prevention		
Handrail	Long handrail (60.63 in and therefore covering the entire belt length)		
Weight and Dimensions (treadmill)		Weight and Dimensions (treadmill incl. packaging)	
Length	82.3 in +/- 1/2 in	Length	90.55 in
Width	33.8 in +/- 1/2 in	Width	43.31 in
Height	51.6 in +/- 1/2 in	Height	38.58 in
Net. weight treadmill	approx. 465 lbs	Gross weight	approx.: 606 – 661 lbs (depending on packaging and transportation)
Subject to technical changes without prior notice.			